

waterlase*MDX™ All-Tissue Laser

**THE BEST-SELLING ALL-TISSUE LASER,
MADE EVEN BETTER!**

New!



BIOLASE | **25 YEARS**
THE WORLD LEADER IN DENTAL LASERS

Revolutionizing Surgery in Dentistry and Medicine

NEW ON THE *waterlase**MDX™

Thousands of dentists have added "no shot, no drill, no pain*" WaterLase treatment to their patient offerings with our MD laser system. Once again we've improved this time-tested technology, adding some of the same advanced features found on our industry-leading WaterLase iPlus!

RR POWER PLANT DELIVERS 300mJ AT 8W OR 450mJ AT 9W FOR 50% FASTER CUTTING



300 mJ Model

The most affordable WaterLase — upgradeable to 450 mJ as your laser skills and applications expand.

450mJ Model

The choice for fast cutting on teeth, soft-tissue and bone in a variety of everyday laser procedures.


UPDATED USER INTERFACE



TWO HANDPIECES FOR MAXIMUM VERSATILITY



EXTREMELY LIGHT AND FLEXIBLE TITANIUM FIBER CABLE HELPS ELIMINATE USER FATIGUE



WATERLASE MDX TECHNICAL SPECS

Wavelength: Er;Cr:YSGG, 2780 nm

Power: 0.1 to 8.0 W (300mJ Model), 0.1 to 9.0 W (450mJ Model)

Pulse Repetition Rates: 10-50 Hz Clinician Selectable

Pulse Energy: 0-300mJ (300mJ Model), 0-450mJ (450mJ Model)

Laser Classification: 4

Operating Voltage: 100-230 VAC

Dimensions: 10.5 x 19 x 32 in (27 x 48 x 91 cm)

Weight: 75lbs (34 kg)

Specifications subject to change without notice.

BIOLASE | **25 YEARS**
THE WORLD LEADER IN DENTAL LASERS

BIOLASE Technology, Inc.

4 Cromwell
Irvine, CA 92618
949.361.1200

Toll-free 888.424.6527
biolase.com
NASDAQ: BLTI

*Discomfort is based on individual sensitivity to pain, and other factors, so not all patients can be treated painlessly without anesthetic.

5201329 Rev. A

