

VistaCam iX HD

the innovative multi-head-camera



مزایای دستگاه Vistacam ix HD :

- . سر قابل تعویض پروکسیمالی
- . انتقال اطلاعات بصورت کاملا دیجیتال
- . استفاده آسان به دلیل شکل ارگونومیک دستگاه
- . قابلیت ضبط تصاویر به صورت عکس یا ویدئو

. کیفیت بالای عکسهای داخل دهانی High Defenition

Technical Data	VistaCam iX HD	VistaCam iX PC version
Connections	USB 2.0 (USB 3.0 compatible)	USB 2.0
Multuser application	Plug & Play	Plug & Play, cable can be connected to handpiece
Activation	Via switch on hand piece (top or bottom, as desired)	Handpiece (optional foot switch)
Handpiece weight	70 g	50 g (without cable)
Handpiece length	200 mm	190mm
Cable length	(optional 3m extension via active holder with USB hub)	2.5 m
Power supply	USB (5V)	USB (5V)
Sensor	High performance CMOS sensor	High performance CMOS sensor
Driver	Uses standard Windows drivers, NO additional drivers needed	Uses standard Windows drivers, NO additional drivers needed
Resolution	1280 pixel (H) x 1024 pixel (V)	
Lighting	2 LEDs, Proof (405 nm, violet), Proxi (850 nm, infrared)	8 LEDs, Cam (whait)4 LEDs, Proof (violet) 2 LEDs,infrared
Optical system	Several lenses and protective glass	Several lenses and protective glass



شرکت بازرگانی نیشام آریانا

www.nisham-aryana.com

تلفن: ۰۲۱-۷۵۳۸۱



high_tech_dent86



nisham_aryana